# Spengler

## TENSONIC®

CE  
0459

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**CHRN-AR-20002390**  
TAS SAT AG  
Chamerstrasse 172, 6300 Zug, Switzerland

## CONTENTS

Safety Notice .....	41
Unit Illustration.....	47
Important Testing Guideline.....	49
Quick Start.....	50
Unit Operation.....	51
Battery Installation.....	51
System Settings.....	52
Applying The Wrist Monitor .....	54
Testing.....	57
Power Off.....	58
Memory Check and test Average.....	59
Memory Deletion.....	60
Low Battery Indicator .....	60
Troubleshooting.....	61
Blood Pressure Information .....	62
Blood Pressure Q&A.....	64
Maintenance.....	65
Specifications.....	66
Warranty.....	69
Electromagnetic Compatibility Information.....	69

## SAFETY NOTICE

Thank you for purchasing the Tensonic® Blood Pressure Monitor. The unit has been constructed using reliable circuitry and durable materials. Used properly, this unit will provide years of satisfactory use.

This device is intended for non-invasive measuring of an adult individuals' systolic, diastolic blood pressure and heart rate using the oscillometric method. The device is not intended for use on infants and children under the age of 12. The device is designed for home or clinical use. All values can be read out in one LCD DISPLAY. Measurement position is on adult wrist only.

Please read this manual thoroughly before using the unit. Please retain this manual for future reference. For specific information about your blood pressure, please CONSULT YOUR DOCTOR. The PATIENT is an intended OPERATOR.

To avoid risk and damage follow all warning precautions. Operate unit only as intended. Read all instructions prior to use.

### MATERIOVIGILANCE

All serious incidents related to the use of a medical device should be notified to the manufacturer and the responsible authority of the Member State where the user and/or the patient live.

## WARNING SIGNS AND SYMBOLS USED

	Caution
	Mandatory
	Prohibited
	Type BF Equipment
	Instructions For Use <b>MUST</b> be Consulted
	Serial Number
	Batch Number
	Discard the used product to the recycling collection point according to local regulations
	The product conforms to the requirements of the EC medical devices directive 93/42/EEC
	Manufacturer
	Swiss representative
	Keep Dry
	Keep off Sunlight
	Power supply: batteries
	Medical device

	Do not use if package is damaged
	Recyclable packaging
IP 22	Protected against solid objects larger than 12.5 mm. Protected against falling water up to 15° from the vertical.
	Transportation and storage temperature range
	Humidity range of transportation and storage

 **Caution**

Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.	
Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.	
Make sure the electronic blood pressure doesn't stop the blood circulation.	
Product is designed for its intended use only. Do not misuse in any way.	
Product is not intended for infants or individuals who cannot express their intentions.	
Do not disassemble or attempt to repair.	
Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.	

This device must not be used in presence of electrosurgery.	
This device must not be used during patient transport.	
Do not apply cuff on a wound, risk of more injuries.	

### **Battery Precautions**

Do not mix new and old batteries simultaneously.

Replace batteries when Low Battery Indicator “” appears on screen.

Be sure battery polarity is correct.

Do not mix battery types. Long-life alkaline batteries are recommended.

Remove batteries from device when not in operation for more than 3 months.

Dispose batteries properly; observe local laws and regulations.

## **Important Instructions Before Use**

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. Contact your physician if test results regularly indicate abnormal readings.
3. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
4. Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.
5. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
6. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant


- should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
7. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.
  8. Too frequent measurements can cause injury to the patient due to blood flow interference.
  9. The cuff should not be applied over a wound as this can cause further injury.
  10. DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
  11. The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
  12. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
  13. A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
  14. Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.
  15. Product is designed for its intended use only. Do not misuse in any way.
  16. Product is not intended for infants or individuals who cannot express their intentions.
  17. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.
  18. Do not disassemble the unit or arm cuff. Do not attempt to repair.
  19. Use only the approved arm cuff for this unit. Use of other arm cuffs may result in incorrect measurement results.
  20. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Make sure to store the blood pressure monitor, children, pets and pests are outside of accessible range.
  21. Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, they may cause incorrect readings and

- interference or become interference source to the device.
22. Do not mix new and old batteries simultaneously
  23. Replace batteries when Low Battery Indicator " appears on screen.  
Replace both batteries at the same time.
  24. Do not mix battery types. Long-life alkaline batteries are recommended.
  25. Remove batteries from device when not in operation for more than 3 months.
  26. Do not insert the batteries with their polarities incorrectly aligned.
  27. Dispose batteries properly; observe local laws and regulations.
  28.  Advising operator that Instruction manual/ Booklet must be consulted .
  29. Do not use the device during transport vehicles for influencing measurement accuracy, such as patient transport in an ambulance or helicopter.
  30. Contains small parts that may cause a choking hazard if swallowed by infants.

## UNIT ILLUSTRATION


### Monitor Unit

EN

## UNIT ILLUSTRATION


### Display

## UNIT ILLUSTRATION

### Contents

1. Monitor Unit

2. Owner's Manual

3. Plastic  
Storage Case


4. Batteries  
(x2)

## IMPORTANT TESTING GUIDELINES

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your wrist level with your heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same wrist, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.


## QUICK START

1. Install batteries. (See Figure A)

2. Remove clothing from the wrist area. (See Figure B)

3. Rest for several minutes prior to testing. Wrap cuff around left wrist. (See Figure C)

4. Sit in a comfortable position and place wrist level with heart. (See Figure D)
5. Press "START/STOP" button to start testing. (See Figure E)

Figure D


Figure E

## UNIT OPERATION

### Battery Installation

Slide battery cover off as indicated by arrow.

Install 2 new AAA alkaline batteries according to polarity.

Close battery cover.

Note:

- 1) Replace batteries when Low Battery Indicator “” appears on screen.
- 2) Batteries should be removed from device when not in operation for an extended period of time.

## System settings

With power off, press "Ser" button to actuate system setting. The Memory Group icon flashes.

### 1. Select memory Group

While in the System setting mode you may accumulate test rescuts into 2 different groups. This allows multiple users to save individual test results (up to 60 memories per group). Press "M" button to choose a group setting. The test results will automatically store in each selected group.

### 2. Time/Date setting

Press "SET" button again to set the Time/Date mode. Set the year first by adjusting the "M" button. Press "SET" button again to confirm current month. Continue setting the day, hour and minute in the same way. Every time the "SET" button is pressed, it will lock in your selection and continue in succession (month, day, hour, minute).

### 3. Time Format setting.

Press "SET" button again to set the time format mode. Set the time format by adjusting the "M" button.

EU means European Time US means U.S Time.



### 4. Saved Settings

While in any setting mode, press "START/STOP" button to turn the unit off. All information will be saved.

Note: Unit will automatically save all information and shut off if left it for 3 minutes.

## Applying The Wrist Monitor

Do not apply over clothing. If wearing a long sleeved shirt, be sure to roll sleeve back to forearm. Apply monitor to wrist as illustrated. Tighten cuff firmly as not to wiggle.

Do not stand while measuring. Sit in a comfortable position with your back supported and your feet flat on the floor, without crossing your legs. Position the center of the cuff at the level of the right auricle of the heart, with the arm supported.


## Testing

### 1.PowerOn

Press and hold "START/STOP button until a beep sounds.

The LCD screen will appear for one second as unit performs a quick diagnosis.

A long tone indicates when unit is ready for testing.

Note: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash "↓" until pressure is stabilized.

## 2. Pressurization

Initial pressure is first pumped to 190 mmHg. If the current user's systolic blood pressure is over 190mmHg, the unit will automatically re-inflate to the proper shelf.

Note: Pressurization will gradually subside and ultimately stop when cuff is not properly applied to the wrist. If this occurs, pressing the START/STOP button to turn the unit off.

## 3. Testing

After cuff Inflation, air will slowly subside as indicated by the corresponding cuff pressure value. A flashing "♥" will appear simultaneously on screen signalling heart beat detection.


Note: Remain relaxed during testing. Avoid speaking or moving body parts.

#### 4. Result Display

The screen will display measurements for systolic and diastolic blood pressure. An indicator representing the current measurement will appear next to the corresponding WHO Classification.

Note: Refer to Page 62-63 for detail WHO Blood Pressure Classification Information.

#### Irregular Heartbeat Indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol "❤️" appears on screen along with measurement results.

Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult your physician if the Irregular Heartbeat Symbol "❤️" frequently appears with your test results.

## 5. Deleting/Storing Test Results

User may delete their current test result due to unfavorable testing conditions or for any other reason. To delete the last test result press the "SET" button after result is displayed. If result is not deleted, it will automatically store by date within the previously configured Memory Group.

Note: be sure the appropriate Memory Group selection is made prior to testing.

If the number of tests surpasses the allotted 60 memories per group, the most recent tests will appear first thus eliminating the oldest readings.

## Power off

The "START/STOP" button can be pressed to turn off the unit in any mode. The unit can turn off the power itself about 3 minutes no operation in any mode.

Safety Precaution: If pressure In cuff becomes too extreme while testing, press the « START/STOP » button to turn power off. The cuff pressure will rapidly dissipate once the unit is off.

## Memory Check and Test Average


With power off, you may check past test results and test average using the "M" buttons. First the average value of all saved measurements of this user memory is display.

Press the "M" button again, and the screen will display the average of last seven days test results for 5:00am - 8:59am measurement.

Press the "M" button again, and the screen will display the average of last seven days test results for afternoon 18:00pm - 19:59pm measurement. .

Press the "M" button again, you may check past test results, the most recent test result can be viewed. Upon activating test results, you can press the "M" buttons to scroll through all test results stored in memory.

Note: If there's no results saved, the systolic Blood Pressure, Diastolic Blood Pressure and Pulse Rate will display "-- --".

## Memory Deletion

Memory for a selected group may be deleted while in Memory Check mode. Press and hold the "SET" button for approximately 3 seconds to delete all memory records from the selected group. The monitor will beep indicating successful deletion and then transfer into testing mode. To delete memories from other groups, please refer to "Select Memory Group" on page 53. Then press the "START/STOP" button to turn the unit off.

Note: Memory cannot be recovered once it has been deleted.

## Low Battery Indicator

The unit will broadcast "Low Battery" when battery life is depleting and unable to inflate cuff for testing. The "⊗" appears simultaneously for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.

## Troubleshooting

<b>Problem</b>	<b>Possible Cause</b>	<b>Solution</b>
Blood pressure results are not within typical range	Cuff is too tight or not properly positioned on the wrist	Firmly reposition cuff on wrist making sure no wiggle is present. (See Page 54)
	Inaccurate test results due to body movement or monitor movement	Sit in a relaxed position placing wrist level with heart. Avoid speaking or moving body parts while testing. (See Page 56)
"Err" displayed	Cuff fails to inflate properly	Make sure hose is properly fastened to cuff and monitor unit.
	Improper operation	Read user manual carefully and re-test properly.
	Pressurization is over 300mmHg	Read user manual carefully and re-test properly.


## BLOOD PRESSURE INFORMATION

### Blood Pressure

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg.) Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.

If these measuring numbers become too high, it means the heart is working harder than it should.


Example: fluctuation within a day (male, 35 years old)

### WHO Blood Pressure Classification Indicator

The Tenstonic® is equipped with a classification indicator based on established guidelines from the World Health Organization.


The chart below (color coded on monitor unit) indicates test results.

█ : Blood Pressure Classification Indicator

## Health Reminder

Hypertension is a dangerous disease that can affect the quality of life. It can lead to a lot of problems including heart failure, kidney failure, and cerebral hemorrhaging. By maintaining a healthy lifestyle and visiting your physician on a regular basis, hypertension and relative diseases are much easier to control when diagnosed in the early stages.

**Note:** Do not be alarmed if an abnormal reading occurs.

A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time each day over an extended period of time. Consult your physician if test results remain abnormal.

## **BLOOD PRESSURE Q&A**

**Q:** What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

**A:** Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

**Note:** Abnormal test results may be caused by:

1. Improper cuff placement

Make sure cuff is snug-not too tight or too loose.

Make sure bottom of the cuff is approximately 1-2cm (1/2") above the elbow joint.

2. Improper body position

Make sure to keep your body in an upright position.

3. Feeling anxious or nervous

Take 2-3 deep breaths, wait a few minutes and resume testing.

**Q:** What causes different readings?

**A:** Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.

**Q:** Should I apply the cuff to the left or right arm? What is the difference?

**A:** Either arm can be used when testing, however, when comparing results, the same arm should be used. Testing on your left arm may provide more accurate results as it is located closer to your heart.

**Q:** What is the best time of day for testing?

**A:** Morning time or any time you feel relaxed and stress free.

## MAINTENANCE


1. Avoid dropping, slamming, or throwing the unit.

2. Avoid extreme temperatures. Do not use outdoors.

3. When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent.

4. Cuff Cleaning: Do not soak cuff in water! Apply a small amount of rubbing alcohol to a soft cloth to clean cuff's surface. Use a damp cloth (water-based) to wipe clean. Allow cuff to dry naturally at room temperature.

5. Do not use petrol, thinners or similar solvents.

6. Remove batteries when not in operation for an extended period of time.

7. Do not disassemble product.

8. It is recommended the performance should be checked every 2 years.

9. Expected service life:  
Approximately three years at 10 tests per day.

10. No service and maintenance while it is in use and maintenance only be performed by service personnel. Service and maintenance require parts, repair, technical support will be provided.

## SPECIFICATIONS

Product Description	Wrist-type Fully Automatic Blood Pressure Monitor	
Model	Tensonic®	
Display	LCD Digital Display - Size:49mmx38mm	
Measurement Method	Oscillometric Method	
Measurement Range	Systolic Pressure	60mmHg-280 mmHg
	Diastolic Pressure	30mmHg-200 mmHg
	Pressure	0mmHg - 300mmHg
	Pulse	30 to 180 Beats/Minute

Measurement Accuracy	Pressure	±3mmHg
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	120 Memories with Date and Time	
Function	WHO Classification Indicator	
	Irregular Heartbeat Detection	
	Average measurements	
Power Source	2 Alkaline Batteries Size AAA	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx. 117g (4 oz.) (Excluding Battery)	
Unit Dimensions	Approx. 77 x 64 x 32mm	
Cuff Circumference	Approx. 75 (W)×300(L) mm Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")	
Operating Environment	Temperature	10°C - 40°C (50°F-104°F)
	Humidity	15% - 93% RH
	Pressure	700hPa-1060hPa
Storage Environment	Temperature	-25°C-70°C (-13°F-158°F)
	Humidity	≤93 % HR
Classification	Internal Powered Equipment, Type BF 	
Ingress Protection Rating	IP22	

Specifications are subject to change without notice.

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0459". This blood pressure monitor also complies with mainly following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety EMC standard:

EN 60601-1-2 Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard:

Electromagnetic Disturbances -Requirements And Tests.

Performance standards:

IEC80601-2-30, Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

EN 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

ISO 81060-2, non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type.

### **Correct Disposal of This Product** (Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed with other household waste at the end of its life.

To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly. When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center. Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

## WARRANTY

The Tensonic® Blood Pressure Monitor is guaranteed for 2-years from the date of purchase. If the Tensonic® does not function properly due to defective components or poor workmanship, we will repair or replace it freely. The warranty does not cover damages to your Tensonic® due to improper handling. Please contact local retailer for details.

## ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

**Table 1**

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emission CISPR 11	Group 1, class B.	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Conducted emission CISPR 11	N/A	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

**Table 2**

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ±4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ±4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A	

Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	N/A	
Voltage dips, short interrup- tions and voltage variations on p- ower supply in- put lines IEC 61000-4-11	< 5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5 secretary	N/A	
Power frequen- cy (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m ; 50 Hz or 60Hz	30 A/m ; 50 Hz or 60Hz	Power frequency magnetic fields should be at levels charactertic of a typical location in a typical commercial or hospital envi- ronment.

**Table 3**

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
RF par conduction CEI 61000-4-6	3V for 0.15-80MHz; 6V in ISM and amateur radio bands between 0.15-80MHz	3V for 0.15-80MHz; 6V in ISM and amateur radio bands between 0.15-80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Radiated RF IEC 61000-4-3	385 MHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 2450 MHz, 28 V (m) 5240 MHz, 5500 MHz, 5785 MHz 9 V/m	385 MHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 2450 MHz, 28 V (m) 5240 MHz, 5500 MHz, 5785 MHz 9V/m	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2,7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[ \frac{3.5}{f_1} \right] \sqrt{P}$	800 MHz to 2,7 GHz $d = \left[ \frac{7}{f_1} \right] \sqrt{P}$
0,01	0,12	0,23
0,1	0,38	0,73
1	1,2	2,3
10	3,8	7,3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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**SPENGLER SAS**